


Nthl1 Rabbit Pab For Wb, Elisa - P78549

Item Number: CM0029716



Introduction

Rabbit polyclonal antibody targeting human NTHL1, validated for WB and ELISA. Detects a 34kDa protein, ideal for base excision repair and oxidative DNA damage research.

[Learn More](#)

Parameter	Value
Product Name	NTHL1 Rabbit pAb
Remarks/Alias	FAP3; NTH1; OCTS3; hNTH1; NTHL1
Species	Human
Gene ID (Human)	4913
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-280 of human NTHL1 (NP_002519.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P78549
Protein Weight	34kDa
Shipping	Ice bag