

Ntn5 Rabbit Polyclonal Antibody For Wb And Elisa - Q8Wtr8

Item Number: CM0016633

Introduction

NTN5 Rabbit Polyclonal Antibody validated for WB and ELISA in human, mouse, rat. Targets netrin-5, a secreted protein essential for motor neuron positioning. Detects ~53kDa band. Q8WTR8.

[Learn More](#)

Parameter	Value
Product Name	NTN5 Rabbit pAb
Remarks/Alias	NTN5; netrin-5
Species	Human
GeneID	126147 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-489 of human NTN5 (NP_665806.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q8WTR8
Protein Weight	53kDa
Shipping	Ice bag