

Nup107 Rabbit Mab - P57740

Item Number: CM0002687



Introduction

Nup107 Rabbit Monoclonal Antibody for Western Blot and ELISA applications. Recognizes human Nup107, a 106 kDa nuclear pore complex protein essential for NPC assembly. Ideal for cell biology and nephrogenesis research.

[Learn More](#)

Parameter	Value
Product Name	Nup107 Rabbit mAb
Remarks/Alias	ODG6; NUP84; NPHS11; GAMOS7; Nup107
Species	Human
GeneID (Human)	57122
GeneID	57122
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700-800 of human Nup107 (P57740)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P57740
Protein Weight	106 kDa
Shipping	Ice bag