

Nup133 Rabbit Pab - Q8Wum0

Item Number: CM0026510



Introduction

NUP133 Rabbit Polyclonal Antibody for WB and ELISA. Detects human NUP133, a 129kDa nucleoporin involved in mRNA export and nephrogenesis. Ideal for nuclear pore complex and kidney development research.

[Learn More](#)

Parameter	Value
Product Name	NUP133 Rabbit pAb
Remarks/Alias	GAMOS8; NPHS18; hNUP133; NUP133
Species	Human
GeneID (human)	55746
GeneID	55746
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 730-930 of human NUP133 (NP_060700.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WUM0
Protein Weight	129kDa
Shipping	Ice bag