

Nup153 Rabbit Polyclonal Antibody For Wb And Elisa - P49790

Item Number: CM0029568



Introduction

Rabbit polyclonal antibody targeting human NUP153, a nuclear pore complex protein. Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for studying nucleocytoplasmic transport and viral integration mechanisms.

[Learn More](#)

Parameter	Value
Product Name	NUP153 Rabbit pAb
Remarks/Alias	N153; HNUP153; NUP153
Species	Human
GeneID (Human)	9972
GeneID	9972
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 350-450 of human NUP153 (NP_005115.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49790
Protein Weight	154kDa
Shipping	Ice bag