


Nup160 Rabbit Polyclonal Antibody For Wb, Elisa - Q12769

Item Number: CM0014259



Introduction

NUP160 rabbit polyclonal antibody validated for WB and ELISA. Targets human nuclear pore complex protein Nup160 (162kDa), a key component of the nuclear pore complex involved in mRNA export. Suitable for nucleocytoplasmic transport research.

[Learn More](#)

Parameter	Value
Product Name	NUP160 Rabbit pAb
Remarks/Alias	NPHS19; NUP160
Species	Human
Gene ID (Human)	23279
Gene ID	23279
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 40-150 of human NUP160 (NP_056046.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q12769
Protein Weight	162kDa
Shipping	Ice bag