


Nup214 Rabbit Polyclonal Antibody For Wb, Ip, Elisa - P35658

Item Number: CM0029792



Introduction

NUP214 rabbit polyclonal antibody (SWISS P35658) validated for WB, IP, ELISA. Targets human nuclear pore complex protein Nup214, involved in nucleocytoplasmic transport and adenovirus capsid disassembly.

[Learn More](#)

Parameter	Value
Product Name	NUP214 Rabbit pAb
Remarks/Alias	CAN; CAIN; IIAE9
Species	Human
Gene ID (Human)	8021
Immunogen	Synthetic peptide (aa 1251-1350 of human NUP214, NP_005076.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human
SWISS	P35658
Protein Weight	214kDa
Shipping	Ice bag