


Nur77 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P22736

Item Number: CM0016274



Introduction

High-quality rabbit polyclonal antibody against human NUR77 (NR4A1), validated for WB, IF/ICC, and ELISA. Ideal for nuclear receptor signaling, inflammation, and cancer research.

[Learn More](#)

Parameter	Value
Product Name	NUR77 Rabbit pAb
Remarks/Alias	HMR; N10; TR3; NP10; GFRP1; NAK-1; NGFIB; NUR77
Species	Human
GeneID (Human)	3164
GeneID	3164
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 200-300 of human NUR77 (NP_002126.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	P22736
Protein Weight	64kDa
Shipping	Ice bag