

Nxn Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q6Dkj4

Item Number: CM0024019

Introduction

NXN Rabbit Polyclonal Antibody validated for WB, IF/ICC, ELISA; reacts with Human, Mouse, Rat. Ideal for Wnt signaling pathway and redox regulation studies. Unconjugated, 48kDa target.

[Learn More](#)

Parameter	Value
Product Name	NXN Rabbit pAb
Remarks/Alias	NRX; RRS2; TRG-4; NXN
Species	Human
Gene ID (Human)	64359
Gene ID	64359
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 100-280 of human NXN (NP_001192248.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q6DKJ4
Protein Weight	48kDa
Shipping	Ice bag