

Occludin Rabbit Polyclonal Antibody For Wb And Elisa - Q16625

Item Number: CM0013617



Introduction

Rabbit polyclonal antibody against human Occludin (UniProt Q16625). Suitable for WB and ELISA, cross-reacts with mouse and rat. Ideal for tight junction and cell permeability studies.

[Learn More](#)

Parameter	Value
Product Name	Occludin Rabbit pAb
Remarks/Alias	BLCPMG; PTORCH1; PPP1R115; Occludin
Species	Human
GeneID (Human)	100506658
GeneID	100506658
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 350-450 of human Occludin (NP_002529.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16625
Protein Weight	59kDa
Shipping	Ice bag