

Ogdhl Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Uld0

Item Number: CM0015803

Introduction

OGDHL rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Detects human, mouse, and rat OGDHL, a mitochondrial 2-oxoglutarate dehydrogenase involved in TCA cycle, cell proliferation inhibition, and apoptosis. Ideal for metabolism and cancer research.

[Learn More](#)

Parameter	Value
Product Name	OGDHL Rabbit pAb
Remarks/Alias	YOBELN; OGDHL
Species	Human
GeneID (Human)	55753
GeneID	55753
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-110 of human OGDHL (NP_060715.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9ULD0
Protein Weight	114kDa
Shipping	Ice bag