


Ogt Rabbit Pab For Wb - O15294

Item Number: CM0018886



Introduction

Rabbit polyclonal antibody against OGT (O-GlcNAc transferase), validated for WB, IHC-P, IF/ICC, ELISA. Suitable for human, mouse, rat samples. Ideal for metabolic, epigenetic, and signaling studies.

[Learn More](#)

Parameter	Value
Product Name	OGT Rabbit pAb
Remarks/Alias	OGT1; HRNT1; MRX106; XLID106; HINCUT-1; O-GLCNAC; OGT
Species	Human
GeneID (Human)	8473
GeneID	8473
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 510-699 of human OGT (NP_858058.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O15294
Protein Weight	117kDa
Shipping	Ice bag