

## Oncostatin M Rabbit Pab For Wb, If/Icc, Elisa - P13725

Item Number: CM0026897



### Introduction

Rabbit polyclonal antibody targeting human Oncostatin M (OSM), suitable for WB, IF/ICC, and ELISA. Cross-reacts with mouse and rat. Useful for cytokine signaling, cell proliferation, and inflammation research.

[Learn More](#)

| Parameter        | Value                   |
|------------------|-------------------------|
| Product Name     | Oncostatin M Rabbit pAb |
| Remarks/Alias    | OSM; Oncostatin M       |
| Species          | Human                   |
| GeneID (Human)   | 5008                    |
| GeneID           | 5008                    |
| Immunogen        | Recombinant protein     |
| Source           | Rabbit                  |
| Category         | Polyclonal Antibodies   |
| Application      | WB, IF/ICC, ELISA       |
| Cross Reactivity | Human, Mouse, Rat       |
| SWISS            | P13725                  |
| Protein Weight   | 28kDa                   |
| Shipping         | Ice bag                 |