

Osarf23 Rabbit Pab - Osarf23

Item Number: CM0020178

Introduction

OsARF23 Rabbit pAb is a polyclonal antibody raised in rabbit against recombinant OsARF23 protein (*Oryza sativa*). Validated for WB and ELISA, it specifically detects OsARF23 (~94kDa) in rice samples.

[Learn More](#)

Parameter	Value
Product Name	OsARF23 Rabbit pAb
Remarks/Alias	arf1; ARF23; OsARF1; OsARF23
Species	<i>Oryza sativa</i>
Gene ID	4350610
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 700-800 of <i>Oryza sativa</i> OsARF23. (Q2R3F5)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
Protein Weight	94kDa
Shipping	Ice bag