

Osgi Rabbit Polyclonal Antibody (Pab) - Q0Jq44

Item Number: CM0019043

Introduction

Rabbit polyclonal antibody targeting OSGI (GI) in *Oryza sativa*. Validated for WB and ELISA, with a predicted molecular weight of 127kDa. Ideal for plant biology and agricultural research studies.

[Learn More](#)

Parameter	Value
Product Name	OsGI Rabbit pAb
Remarks/Alias	GI; OsGI
Species	Oryza sativa
Gene ID	4325329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 820-1090 of Oryza sativa OSGI (Q9AWL7)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q0JQ44
Protein Weight	127kDa
Shipping	Ice bag