

Osgsr1 Rabbit Polyclonal Antibody For Wb And Elisa - Q53An3

Item Number: CM0021721

Introduction

Rabbit polyclonal antibody targeting rice OsGSR1 (Q53AN3). Suitable for Western blot (WB) and ELISA. Cross-reacts with Oryza sativa. Immunogen: recombinant protein. Predicted molecular weight: 11 kDa.

[Learn More](#)

Parameter	Value
Product Name	OsGSR1 Rabbit pAb
Remarks/Alias	Gsr1; OJ1001_B06.38; OsGSR1
Species	Oryza sativa
Gene ID	4340712
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of Oryza sativa OsGSR1. (Q53AN3).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q53AN3
Protein Weight	11kDa
Shipping	Ice bag