

OsSERK2 Rabbit Polyclonal Antibody For Wb And Elisa - Q7Xv05

Item Number: CM0023757



Introduction

Rabbit polyclonal antibody targeting *Oryza sativa* OsSERK2, raised against recombinant protein. Validated for WB and ELISA, cross-reacts with rice. Suitable for plant immunity, somatic embryogenesis, and brassinosteroid signaling research.

[Learn More](#)

Parameter	Value
Product Name	OsSERK2 Rabbit pAb
Remarks/Alias	OsSERK1; OsSERK2
Species	<i>Oryza sativa</i>
Gene ID	4336035
Immunogen	Recombinant protein of <i>Oryza sativa</i> OsSERK2
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	<i>Oryza sativa</i>
SWISS	Q7XV05
Protein Weight	70kDa
Shipping	Ice bag