

Oswrky53 Rabbit Pab For Wb, Elisa - Q5W6M7

Item Number: CM0022747

Introduction

OsWRKY53 Rabbit pAb, a polyclonal antibody generated against recombinant Oryza sativa OsWRKY53 (aa 1-181), validated for WB and ELISA. Detects the WRKY transcription factor implicated in abscisic acid and gibberellin signaling, localizing to the nucleus.

[Learn More](#)

Parameter	Value
Product Name	OsWRKY53 Rabbit pAb
Remarks/Alias	OsWRKY53
Species	Oryza sativa
GeneID	4338474
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-181 of Oryza sativa OsWRKY53. (Q5W6M7)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Oryza sativa
SWISS	Q5W6M7
Protein Weight	51kDa
Shipping	Ice bag