


P23/Ptges3 Rabbit Mab For Wb, If/Icc, Elisa - Q15185

Item Number: CM0004960



Introduction

P23/PTGES3 Rabbit Monoclonal Antibody validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Targets prostaglandin E synthase 3, a cytosolic enzyme and chaperone involved in inflammation and transcriptional regulation.

[Learn More](#)

Parameter	Value
Product Name	P23/PTGES3 Rabbit mAb
Remarks/Alias	P23; TEBP; cPGES; P23/PTGES3
Species	Human
GeneID (human)	10728
GeneID	10728
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human P23/PTGES3 (NP_006592.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15185
Protein Weight	19 kDa
Shipping	Ice bag