

P2Rx4 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q99571

Item Number: CM0029228



Introduction

Human P2RX4 rabbit polyclonal antibody, targeting P2X purinoceptor 4. Validated for WB, IF/ICC, ELISA with cross-reactivity to human and mouse. Ideal for ion channel, immunology, and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	P2RX4 Rabbit pAb
Remarks/Alias	P2X4; P2X4R; P2RX4
Species	Human
Gene ID (Human)	5025
Gene ID	5025
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 55-338 of human P2RX4 (NP_002551.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q99571
Protein Weight	43kDa
Shipping	Ice bag