


P38 Mapk Rabbit Pab - Q16539

Item Number: CM0015510



Introduction

p38 MAPK Rabbit pAb is a rabbit polyclonal antibody raised against a synthetic peptide within human p38 MAPK (aa 261-360). Suitable for WB, IHC-P, IF/ICC, and ELISA, with cross-reactivity to human, mouse, and rat. Ideal for MAPK signaling research.

[Learn More](#)

Parameter	Value
Product Name	p38 MAPK Rabbit pAb
Remarks/Alias	RK; p38; CSBP; EXIP; Mxi2; CSBP1; CSBP2; CSPB1; PRKM14; PRKM15; SAPK2A; p38ALPHA; p38 MAPK
Species	Human
GeneID (Human)	1432
GeneID	1432
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 261-360 of human p38 MAPK (NP_620581.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q16539
Protein Weight	41kDa
Shipping	Ice bag