


P70 S6 Kinase 1 Rabbit Polyclonal Antibody For Wb - P23443

Item Number: CM0020316



Introduction

Rabbit polyclonal antibody targeting human p70 S6 Kinase 1 (S6K1), validated for WB, IHC-P, IF/ICC, ELISA. Suitable for cell growth and translation regulation research in human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	p70 S6 Kinase 1 Rabbit pAb
Remarks/Alias	S6K; P56K; S6K1; STK14A; p70-S6K; p70 S6KA; p70-alpha; S6K-beta-1; p70(S6K)-alpha; p70 S6 Kinase 1
Species	Human
GeneID	6198
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 426-525 of human P70 S6K (NP_003152.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P23443
Protein Weight	59kDa
Shipping	Ice bag