


P75Ntr/Ngfr Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P08138

Item Number: CM0016814



Introduction

Rabbit polyclonal antibody targeting human p75NTR/NGFR (CD271; TNFRSF16). Validated for WB, IF/ICC, ELISA; cross-reacts with mouse and rat. Suitable for studying neuronal differentiation, survival, and insulin-dependent glucose uptake.

[Learn More](#)

Parameter	Value
Product Name	p75NTR/NGFR Rabbit pAb
Remarks/Alias	CD271; p75NTR; TNFRSF16; p75(NTR); Gp80-LNGFR; p75NTR/NGFR
Species	Human
Gene ID (Human)	4804
Gene ID	4804
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 273-427 of human p75NTR/NGFR (NP_002498.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08138
Protein Weight	45kDa
Shipping	Ice bag