


P90Rsk/Rsk1/Rps6Ka1 Rabbit Pab - Q15418

Item Number: CM0024648



Introduction

Rabbit polyclonal antibody against p90Rsk/RSK1/RPS6KA1, validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. 83kDa protein. Ideal for MAPK signaling and cell proliferation studies.

[Learn More](#)

Parameter	Value
Product Name	p90Rsk/RSK1/RPS6KA1 Rabbit pAb
Remarks/Alias	RSK; HU-1; RSK1; p90Rsk; MAPKAPK1; MAPKAPK1A; p90Rsk/RSK1/RPS6KA1
Species	Human
GeneID	6195
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 600 to the C-terminus of human p90Rsk/RSK1/RPS6KA1 (NP_002944.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q15418
Protein Weight	83kDa
Shipping	Ice bag