

## Pai-1/Serpin E1 Rabbit Mab For Wb, Elisa - P05121

Item Number: CM0002309



### Introduction

PAI-1/Serpin E1 Rabbit Monoclonal Antibody validated for WB and ELISA. Reacts with human and monkey PAI-1. Ideal for fibrinolysis, cell migration, and senescence research.

[Learn More](#)

| Parameter        | Value                                                          |
|------------------|----------------------------------------------------------------|
| Product Name     | PAI-1/Serpin E1 Rabbit mAb                                     |
| Remarks/Alias    | PAI; PAI1; PAI-1; PLANH1; PAI-1/Serpin E1                      |
| Species          | Human                                                          |
| GeneID (Human)   | 5054                                                           |
| Immunogen        | Synthetic Peptide (aa 50-150 of human PAI-1/Serpin E1, P05121) |
| Source           | Rabbit                                                         |
| Category         | Monoclonal Antibodies                                          |
| Application      | WB, ELISA                                                      |
| Cross Reactivity | Human, Monkey                                                  |
| SWISS            | P05121                                                         |
| Protein Weight   | 45kDa                                                          |
| Shipping         | Ice bag                                                        |