

Pak1 Rabbit Mab For Wb, If/Icc, Elisa - Q13153

Item Number: CM0002644



Introduction

Rabbit monoclonal antibody targeting human PAK1 (Q13153), validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for cancer, cytoskeleton, and signaling pathway studies.

[Learn More](#)

Parameter	Value
Product Name	PAK1 Rabbit mAb
Remarks/Alias	IDDMSSD; p65-PAK; PAKalpha; alpha-PAK; PAK1
Species	Human
GeneID (human)	5058
GeneID	5058
Immunogen	Synthetic Peptide (A synthetic peptide corresponding to a sequence within amino acids 1-100 of human PAK1 (Q13153))
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13153
Protein Weight	61kDa
Shipping	Ice bag