

Pak1 Rabbit Polyclonal Antibody For Wb And Elisa - Q13153

Item Number: CM0015009



Introduction

High-quality rabbit polyclonal antibody targeting human PAK1 (alpha-PAK, p65-PAK). Validated for WB and ELISA, cross-reacts with rat. Ideal for studying cytoskeleton dynamics, cell migration, and apoptosis. Recombinant immunogen ensures specificity.

[Learn More](#)

Parameter	Value
Product Name	PAK1 Rabbit pAb
Remarks/Alias	IDDMSSD; p65-PAK; PAKalpha; alpha-PAK; PAK1
Species	Human
GeneID (Human)	5058
GeneID	5058
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 130-250 of human PAK1 (NP_001122092.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	Q13153
Protein Weight	61kDa
Shipping	Ice bag