


Pak5 Rabbit Polyclonal Antibody For Wb, Elisa - Q9P286

Item Number: CM0022242



Introduction

PAK5 Rabbit Polyclonal Antibody validated for Western Blotting (WB) and ELISA, with reactivity in human and mouse. Ideal for studying PAK5/PAK7 in cytoskeletal regulation, cell signaling, and brain research.

[Learn More](#)

Parameter	Value
Product Name	PAK5 Rabbit pAb
Remarks/Alias	PAK7; PAK5
Species	Human
GeneID (human)	57144
GeneID	57144
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 600-700 of human PAK5 (NP_065074.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9P286
Protein Weight	81kDa
Shipping	Ice bag