

Pan-Akt Rabbit Mab For Western Blot, Ip, If, Ihc-P, Elisa - P31749 / P31751 / Q9Y243

Item Number: CM0002895



Introduction

Pan-Akt Rabbit Monoclonal Antibody validated for WB, IP, IF, IHC, ELISA. Detects Akt1, Akt2, Akt3 in human, mouse, rat. Molecular weight: 56kDa.

[Learn More](#)

Parameter	Value
Product Name	Pan-Akt Rabbit mAb
Remarks/Alias	AKT; PKB; RAC; PRKBA; PKB-ALPHA; RAC-ALPHA
Species	Human
GeneID (Human)	207, 208, 10000
GeneID	207, 208, 10000
Immunogen	Recombinant protein: Recombinant fusion protein containing a sequence corresponding to amino acids 1-123 of human Pan-Akt (NP_005154.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749 / P31751 / Q9Y243
Protein Weight	56kDa
Shipping	Ice bag