

Pan-Akt Rabbit Polyclonal Antibody For Wb, Ip, If/Icc, Ihc-P, Elisa - P31749/P31751/Q9Y243

Item Number: CM0022311



Introduction

Rabbit polyclonal antibody targeting pan-Akt (AKT1, AKT2, AKT3) for WB, IP, IF/ICC, IHC-P, and ELISA. Cross-reacts with human, mouse, rat. Ideal for studying AKT pathways in cancer and metabolism.

[Learn More](#)

Parameter	Specification
Product Name	Pan-Akt Rabbit pAb
Remarks/Alias	AKT1/AKT2/AKT3; Pan-Akt
Species	Human
GeneID (Human)	207/208/10000
GeneID	207/208/10000
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 311-480 of human AKT1/AKT2/AKT3 (NP_005154.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31749/P31751/Q9Y243
Protein Weight	48kDa/55kDa/51kDa/54kDa
Shipping	Ice bag