

Pan-Crotonyl-K Rabbit Mab For Wb, If/Icc, Elisa - High Specificity Lysine Crotonylation Antibody

Item Number: CM0000017

Introduction

Pan-Crotonyl-K Rabbit monoclonal antibody (mAb) validated for WB, IF/ICC, ELISA. Recognizes crotonylated lysine residues across human, mouse, rat, and other species. Ideal for studying protein crotonylation modifications.

[Learn More](#)

Parameter	Value
Product Name	Pan-Crotonyl-K Rabbit mAb
Immunogen	Synthetic peptide
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Protein Weight	10kDa or more
Shipping	Ice bag