

Pan-Crotonyl-K Rabbit Polyclonal Antibody For Wb And Elisa

Item Number: CM0000030

Introduction

Rabbit polyclonal antibody specific for crotonylated lysine residues. Validated for Western blot and ELISA, reacts with human, mouse, and rat samples. Ideal for studying lysine crotonylation modifications.

[Learn More](#)

Parameter	Details
Product Name	Pan-Crotonyl-K Rabbit pAb
Species	Human
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence containing crotonylated K.
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag