

Pan Dimethyl-Lysine Rabbit Pab For Wb, Ihc-P, If/Icc, Elisa

Item Number: CM0000227

Introduction

Pan DiMethyl-lysine Rabbit Polyclonal Antibody (pAb) targeting dimethylated lysine residues. Suitable for Western blot, IHC-paraffin, IF/ICC, and ELISA. Predicted to cross-react with human, mouse, rat, and other species. Immunogen: synthetic peptide containing dimethylated K.

[Learn More](#)

Parameter	Value
Product Name	Pan DiMethyl-lysine Rabbit pAb
Gene ID (Human)	pan Kme2
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence containing dimethylated K.
Source	Rabbit
Category	Methylated Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag