

Pan Lactic Acid-Lysine Rabbit Monoclonal Antibody For Wb, Ip, Elisa | Lactylation Detection

Item Number: CM0000061

Introduction

Rabbit monoclonal antibody targeting Pan Lactic acid-Lysine, a key post-translational modification. Validated for WB, IP, ELISA with broad species reactivity (Human, Mouse, Rat). Ideal for lactylation research.

[Learn More](#)

Parameter	Value
Product Name	Pan Lactic acid-Lysine Rabbit mAb
Species	Human
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence containing Pan Lactic acid-Lysine.
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag