

Pan Lactic Acid-Lysine Rabbit Polyclonal Antibody (Wb, Elisa)

Item Number: CM0000053

Introduction

Rabbit polyclonal antibody targeting Pan Lactic acid-Lysine modification. Validated for Western blot (WB) and ELISA applications. Reacts with human, mouse, rat, and other wide-range species. Ideal for studying protein post-translational modifications.

[Learn More](#)

Parameter	Value
Product Name	Pan Lactic acid-Lysine Rabbit pAb
Species	Human
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence containing Pan Lactic acid-Lysine
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range)
Shipping	Ice bag