

Pan-Mono-Methyl Arginine Motif Rabbit Polyclonal Antibody For Wb And Elisa

Item Number: CM0000224

Introduction

Rabbit polyclonal antibody targeting pan-mono-methyl arginine motifs. Validated for Western blot (WB) and ELISA applications, cross-reactive with human, mouse, and rat. Suitable for arginine methylation studies.

[Learn More](#)

Parameter	Value
Product Name	pan-Mono-Methyl Arginine Motif Rabbit pAb
GeneID (Human)	pan Rme1
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence containing monomethylated R.
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag