

Pan Mono-Methyl Lysine Rabbit Polyclonal Antibody For Western Blot And Elisa Applications

Item Number: CM0000213

Introduction

Pan Mono-Methyl Lysine Rabbit Polyclonal Antibody (Cat. No. A18293) for Western Blot and ELISA. Reacts with human, mouse, and rat. Ideal for protein methylation research. Reliable sourcing from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	Pan Mono-Methyl lysine Rabbit pAb
Gene ID (Human)	pan Kme1
Immunogen	`Synthetic peptide
Source	Rabbit
Category	Methylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag