

Pan-Mono-Methylr (R*Gg) Motif Rabbit Polyclonal Antibody For Wb, Db, Elisa

Item Number: CM0000219

Introduction

Validated rabbit polyclonal antibody targeting pan-mono-methylarginine R*GG motif for Western blot, dot blot, and ELISA. Reacts with human, mouse, rat. Ideal for global methylation detection research.

[Learn More](#)

Parameter	Value
Product Name	pan-Mono-MethylR R*GG Motif Rabbit pAb
Gene ID (Human)	pan Rme1
Immunogen	Synthetic peptide
Source	Rabbit
Category	Methylated Antibodies
Application	WB, DB, ELISA
Cross Reactivity	Human, Mouse, Rat
Shipping	Ice bag