

Pan-Myristyl-G Rabbit Polyclonal Antibody For Wb And Elisa

Item Number: CM0000006

Introduction

pan-Myristyl-G Rabbit Polyclonal Antibody (pAb) targets myristylated glycine modifications. Validated for Western Blot (WB) and ELISA applications, with rat cross-reactivity. Ideal for studying protein lipidation and myristoylation in research.

[Learn More](#)

Parameter	Value
Product Name	pan-Myristyl-G Rabbit pAb
Target	pan Myristyl-G
Immunogen	Synthetic peptide; A synthetic peptide corresponding to a sequence containing myristylated G.
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
Shipping	Ice bag