

Pan-Nitrific-Y Rabbit Polyclonal Antibody For Western Blot And Elisa

Item Number: CM0000036

Introduction

Rabbit polyclonal antibody targeting pan-Nitrific-Y (tyrosine nitrication) for Western blot (WB) and ELISA applications. Predicted wide species cross-reactivity with human, mouse, rat, and others. Immunogen: synthetic peptide containing myristylated lysine.

[Learn More](#)

Parameter	Value
Product Name	pan-Nitrific-Y Rabbit pAb
Species	Human
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence containing myristylated K.
Source	Rabbit
Category	Other Modified Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag