

Pan Phospho-Serine/Threonine Mouse Mab - Pan Phospho-Serine/Threonine

Item Number: CM0001561

Introduction

This pan phospho-Serine/Threonine mouse monoclonal antibody is designed for detecting phosphorylated proteins in various species including human, mouse, and rat. It is suitable for western blotting, immunohistochemistry, and ELISA applications in signaling research.

[Learn More](#)

Parameter	Value
Product Name	pan Phospho-Serine/Threonine Mouse mAb
Species	Universal
Gene ID (Human)	pan p-S/T
Immunogen	A synthetic peptide corresponding to a sequence containing phosphorylated Serine/Threonine.
Source	Mouse
Category	Phosphorylated Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag