

Pan Phospho-Serine/Threonine Rabbit Mab - Phospho-Serine/Threonine Antibody

Item Number: CM0000844

Introduction

Pan-specific rabbit monoclonal antibody targeting phosphorylated serine/threonine residues. Suitable for WB, IP, IHC-P, and ELISA across human, mouse, rat, and other predicted species. Ideal for broad phospho-signaling studies.

[Learn More](#)

Parameter	Value
Product Name	pan Phospho-Serine/Threonine Rabbit mAb
Species	Universal
Immunogen	Synthetic peptide corresponding to a sequence containing phosphorylated Serine/Threonine
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag