

Pan Phospho-Serine/Threonine Rabbit Polyclonal Antibody - Pan P-S/T

Item Number: CM0000545

Introduction

Pan Phospho-Serine/Threonine Rabbit Polyclonal Antibody for detection of phosphorylated S/T residues. Validated for WB and ELISA applications, with predicted wide cross-reactivity across human, mouse, rat, and other species. Ideal for studying phosphorylation-dependent signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	Pan Phospho-Serine/Threonine Rabbit pAb
Gene ID	pan p-S/T
Immunogen	A synthetic peptide corresponding to a sequence containing phosphorylated S & T.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag