

Pan Phospho-T Rabbit Pab - Phosphothreonine (Phospho-T)

Item Number: CM0000675

Introduction

Rabbit polyclonal antibody targeting phosphorylated threonine residues, suitable for Western blot and ELISA. Predicted to react with a wide range of species including human and mouse, useful for detecting protein phosphorylation events in various research contexts.

[Learn More](#)

Parameter	Value
Product Name	Pan Phospho-T Rabbit pAb
Immunogen	Synthetic peptide corresponding to a sequence containing phosphorylated T
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Other (Wide Range Predicted)
Shipping	Ice bag