

Pan Phospho-Threonine Rabbit Mab - Pan Phospho-Threonine

Item Number: CM0001394

Introduction

Pan Phospho-Threonine Rabbit mAb detects threonine phosphorylation across a wide range of species (human, mouse, rat, and more). This rabbit monoclonal antibody is validated for WB, IP, and ELISA, supporting phosphorylation signaling research.

[Learn More](#)

Parameter	Value
Product Name	Pan Phospho-Threonine Rabbit mAb
Species	Human
Immunogen	A synthetic phosphorylated peptide around of human Threonine.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag