

Pan Phospho-Tyrosine Mouse Mab - Pan P-Y

Item Number: CM0001112

Introduction

Pan Phospho-Tyrosine Mouse monoclonal antibody detects phosphorylated tyrosine residues. Suitable for WB, IF/ICC, ELISA applications with predicted reactivity across human, mouse, rat, and a wide range of species. Essential for tyrosine kinase signaling research.

[Learn More](#)

Parameter	Value
Product Name	Pan Phospho-Tyrosine Mouse mAb
Gene ID (Human)	pan p-Y
Immunogen	A synthetic peptide corresponding to a sequence containing phosphorylated Y
Source	Mouse
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag