

Pan Phospho-Tyrosine Rabbit Mab - Pan P-Y

Item Number: CM0001604

Introduction

Rabbit monoclonal antibody recognizing pan-phospho-tyrosine residues across human, mouse, and rat species, validated for WB, IP, and ELISA applications. Ideal for studying tyrosine phosphorylation signaling networks.

[Learn More](#)

Parameter	Value
Product Name	Pan Phospho-Tyrosine Rabbit mAb
Species	Human
Gene ID (Human)	Pan p-Y
Immunogen	Synthetic Peptide: A random peptide library phosphorylated peptide around Tyr
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag