

Pan Phospho-Tyrosine Rabbit Monoclonal Antibody (Mab) - Pan Phospho-Tyrosine Rabbit Mab

Item Number: CM0001976

Introduction

Pan Phospho-Tyrosine Rabbit monoclonal antibody (mAb) for WB, IP, ELISA. Detects phospho-tyrosine across human, mouse, rat, and other species with broad reactivity. Ideal for phosphorylation and cell signaling studies. >10KDa.

[Learn More](#)

Parameter	Value
Product Name	Pan Phospho-Tyrosine Rabbit mAb
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence containing phosphorylated Y.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IP, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Protein Weight	>10KDa
Shipping	Ice bag