

Pan Phospho-Tyrosine Rabbit Pab - Pan P-Y

Item Number: CM0001258

Introduction

Rabbit polyclonal antibody targeting pan phospho-tyrosine, suitable for WB, IF/ICC, ELISA applications. Recognizes human, mouse, rat, and a wide range of other species. Ideal for tyrosine phosphorylation research.

[Learn More](#)

Parameter	Value
Product Name	Pan Phospho-Tyrosine Rabbit pAb
Gene ID (Human)	Pan p-Y
Immunogen	Synthetic Peptide. A synthetic peptide corresponding to a sequence containing phosphorylated S & T.
Source	Rabbit
Category	Phosphorylated Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat, Other (Wide Range Predicted)
Shipping	Ice bag