

## Pan-Succinyl-K Rabbit Mab - Succinylated Lysine Antibody For Wb, Ip, Db, Elisa

Item Number: CM0000013



### Introduction

Pan-Succinyl-K Rabbit monoclonal antibody for detecting succinylated lysine residues. Suitable for WB, IP, DB, ELISA with predicted wide cross-reactivity (human, mouse, rat). Ideal for post-translational modification research.

[Learn More](#)

| Parameter        | Value                                                                                         |
|------------------|-----------------------------------------------------------------------------------------------|
| Product Name     | Pan-Succinyl-K Rabbit mAb                                                                     |
| Immunogen        | Synthetic peptide: A synthetic peptide corresponding to a sequence containing succinylated K. |
| Source           | Rabbit                                                                                        |
| Category         | Other Modified Antibodies                                                                     |
| Application      | WB, IP, DB, ELISA                                                                             |
| Cross Reactivity | Human, Mouse, Rat, Other (Wide Range Predicted)                                               |
| Shipping         | Ice bag                                                                                       |